

HYDROGEL 23 - PRODUCTSHEET



HYDROGEL 23 – UPPERLEG COOLING BRACE

The IZI hydrogel 23 upperleg cooling brace combines function, comfort, and a great fit to reduce pain at the upperleg. Our upperleg braces are lightweight, low-profile, with soft neoprene, and easily adjustable by its velcro straps . Our self-cooling upperleg brace includes our patented hydrogel cooling technologie . IZI cooling braces are unique in their kind due to their double function the top choice for upperleg cooling , support and relief from pain at the upperleg.

The upperleg cooling brace is activated by putting it in a bucket of (cold) water for 20 minutes . For additional cooling , the upperleg brace can be frozen. This gives about half hour of intense cooling , which is comparable to cooling with ice packs. When the ice has melted, the upperleg brace will continue the cooling by evaporation .

MATERIALS USED :

FABRIC 1 WOVEN POLYESTER, 85 GRMS. (COOLING PART)
 FABRIC 2 NEOPREEN (blue color)

FINISHINGS FABRIC 1 + QUICK DRY FINISH
 FABRIC 1 + CIRE FINISH

AVAILABLE COLORS BASED ON ABOVE FINISHINGS :



NOTE : THE PRODUCT ON THIS SHEET IS REFERENCE FOR COLORS IN STOCK.

OTHER COOLING BRACES IZI HAS DEVELOPED :

- | | | |
|-------------|---|--------------------|
| HYDROGEL 20 | > | WRIST (LEFT HAND) |
| HYDROGEL 21 | > | WRIST (RIGHT HAND) |
| HYDROGEL 22 | > | ELBOW |
| HYDROGEL 24 | > | LOWERLEG |
| HYDROGEL 25 | > | ANKLE |
| HYDROGEL 26 | > | KNEE |

FEATURES :

WEIGHT 50 GRMS. (DRY)
 150 GRMS. (ACTIVATED)
 COOLING TIME APPROX. 72 HOURS EVAPORATED COOLING
 APPROX. 2 HOURS AS ICE BRACE
 COOLING METHOD EVAPORATION AND ICE
 FITTING ADJUSTABLE
 SIZES 1 SIZE FITS ALL
 MEDICAL CE YES

Printed size indication



**1 SIZE
 FITS ALL**

HYDROGEL COOLPANEL BACK



SOFT NEOPRENE FOR COMFORT



VELCRO FOR ADJUSTABLE
 FIT AT WRIST