

HYDROGEL 26 - PRODUCTSHEET



HYDROGEL 26 – KNEE COOLING BRACE

The IZI hydrogel 26 knee cooling brace combines function, comfort, and a great fit to reduce pain at the knee. Our knee braces are lightweight, low-profile, with soft neoprene, and easily adjustable by its velcro straps . Our self-cooling knee brace includes our patented hydrogel cooling technologie . IZI cooling braces are unique in their kind due to their double function the top choice for knee cooling , support and relief from pain at the knee.

The knee cooling brace is activated by putting it in a bucket of (cold) water for 20 minutes . For additional cooling , the knee brace can be frozen. This gives about half hour of intense cooling , which is comparable to cooling with ice packs. When the ice has melted, the knee brace will continue the cooling by evaporation .

MATERIALS USED :

FABRIC 1 WOVEN POLYESTER, 85 GRMS. (COOLING PART)
FABRIC 2 NEOPREEN (blue color)

FINISHINGS FABRIC 1 + QUICK DRY FINISH
FABRIC 1 + CIRE FINISH

AVAILABLE COLORS BASED ON ABOVE FINISHINGS :



NOTE : THE PRODUCT ON THIS SHEET IS REFERENCE FOR COLORS IN STOCK.

OTHER COOLING BRACES IZI HAS DEVELOPED :

| | | |
|-------------|---|--------------------|
| HYDROGEL 20 | > | WRIST (LEFT HAND) |
| HYDROGEL 21 | > | WRIST (RIGHT HAND) |
| HYDROGEL 22 | > | ELBOW |
| HYDROGEL 23 | > | UPPERLEG |
| HYDROGEL 24 | > | LOWERLEG |
| HYDROGEL 25 | > | ANKLE |

FEATURES :

| | |
|----------------|---|
| WEIGHT | 50 GRMS. (DRY) 150 GRMS. (ACTIVATED) |
| COOLING TIME | APPROX. 72 HOURS EVAPORATED COOLING APPROX. 2 HOURS AS ICE BRACE |
| COOLING METHOD | EVAPORATION AND ICE |
| FITTING | ADJUSTABLE |
| SIZES | 1 SIZE FITS ALL |
| MEDICAL CE | YES |

Printed size indication



**1 SIZE
FITS ALL**

HYDROGEL COOLPANEL BACK



SOFT NEOPRENE FOR COMFORT



VELCRO FOR ADJUSTABLE
FIT AT WRIST



HYDROGEL
IZI COOLING TECHNOLOGY